

MediEase Technical Requirements

Network Server

Operating System

- MediEase runs on Microsoft NT and all Novell Networks. (Novell 3.51 would require adjusting long file names).
- All users require read, write and create user rights to the MediEase installation and data directories.
- Ability to access SIS export file location

Storage

- 10 gigs should provide storage for up to 3.5 million patient records. A safe rule of thumb is 25 megabytes for the initial installation with all programs, help and utility files and 30 megabytes for each 10,000 patient records.

Work Stations

General

- Microsoft Windows 98, 2000 or NT operating systems. Processor should be at least a 486. Pentium is recommended (especially for reporting) with the minimum speed of 300mhz
- A minimum of 64 megabytes of RAM (96 MB or more are recommended)
- Compatible VGA or higher resolution display supported by Microsoft Windows
- Mouse
- 32 megabytes minimum of free disk space on the C: drive

For Reporting

- 500 megabytes of free disk space on the C: drive.

Administrator

- Internet access suggested for downloading updates

Other

- Printers should be compatible with Microsoft Windows - no special requirements
- If using email to send letters, the user will need access to the Internet.